

The auditory hemodynamic response & Loudness modulation

Parzijn Majid 1:50 AM

good afternoon teammates, I am a clinical linguist and founder of earlab co. Earlab is a start-up company. we want to develop a medical diagnostic instrument for measuring hearing test. the instrument is based on measuring brain activity using EEG in combination with fnirs. this has also been the reason for me to join this group.

I don't know much about BCI or programming and no experience with fnirs jet!.

Parzijn Majid 1:56 AM

We are in the early Development phase(POC)

One of our assumptions that we are working on is Loudness modulation for different frequency bands (between 100 Hz to 20kHz) in an interval of a third (1/3 octave)

Lawrence Neal 3:42 AM

FYI, I'm a Retired Telecommunications Engineer - 71yo, active as an Embedded Systems Developer, Hardware/Software, currently working on a project titled 'EpiSense' to anticipate Aura/Seizures for Epilepsy, A. Generalized tonic-clonic, B. Dyscognitive, C. Generalized tonic-clonic, and D. Dyscognitive seizures.

Yes, we need "Hardware", preferably the Nautilus fNIRS at minimum.

Parzijn Majid 4:00 AM

Most studies performed with fMRI have demonstrated that the auditory hemodynamic response is sensitive to variations in sound level. Some authors have indicated a positive, nearly linear relationship between the strength of the BOLD signal and sound intensity. It appears that auditory cortical responses measured with fNIRS show such a linear relationship for perceived loudness, rather than for the (physical) intensity of the sound. This potential discrepancy between intensity vs. loudness might suggest that fNIRS does not primarily target primary auditory cortex, where intensity effects seem more clear, but mainly relate to activity generated in secondary auditory areas

Greeshma Sharma 4:11 AM

So we are developing BCI application for auditory stimulation.

What protocol should we look after

<https://www.frontiersin.org/articles/10.3389/fnins.2014.00373/full>

FrontiersFrontiers

Hybrid fNIRS-EEG based classification of auditory and visual perception processes

For multimodal Human-Computer Interaction (HCI), it is very useful to identify the modalities on which the user is currently processing information. This would enable a system to select complementary output modalities to reduce the user's workload. In this paper, we develop a hybrid Brain-Computer Interface (BCI) which uses Electroencephalography (EEG) and functional Near Infrared Spectroscopy (fNIRS) to discriminate and detect visual and auditory stimulus processing. We describe the experimental setup we used for collection of our data corpus with 12 subjects. On this data, we performed cross-validation evaluation, of which we report accuracy for different classification conditions. The results show that the subject-dependent systems achieved a classification accuracy of 97.8% for discriminating visual and auditory perception processes from each other and a classification accuracy of up to 94.8% for detecting modality-specific processes independently of other cognitive activity. The same classification conditions could also be discriminated in a subject-independent fashion with accuracy of up to 94.6 and 86.7%, respectively. We also look at the contributions of the two signal types and show that the fusion of classifiers using different features significantly increases accuracy.

Parzijn Majid 4:15 AM

<https://www.degruyter.com/document/doi/10.1515/cdbme-2016-0051/html>

De GruyterDe Gruyter

fNIRS for future use in auditory diagnostics

Functional near-infrared spectroscopy (fNIRS) is an emerging technique for the assessment of functional activity of the cerebral cortex. Recently fNIRS was also envisaged as a novel neuroimaging approach for measuring the auditory cortex (AC) activity in cochlear implant (CI) users. In the present study we report on initial measurements of AC activation due to spatial sound presentation with a first target to generate data for comparison with CI user and the future use in auditory diagnostics.

Parzijn Majid 6:02 AM

<https://www.ncbi.nlm.nih.gov/pmc/articles/PMC6394942/>

PubMed Central (PMC)PubMed Central (PMC)

Assessing hearing by measuring heartbeat: The effect of sound level

Functional near-infrared spectroscopy (fNIRS) is a non-invasive brain imaging technique that measures changes in oxygenated and de-oxygenated hemoglobin concentration and can provide a measure of brain activity. In addition to neural activity, fNIRS signals ...

Default Setting



fNIRS
Project

OFF TIME rest	30 sec	30 sec	30 sec
random off time	Do, Re, Mi,..... 15 sec	Do, Re, Mi,..... 15 sec	Do, Re, Mi,..... 15 sec
rest	10 ms	10 ms	10 ms
tone length signal	5 sec Pure tone	5 sec Pure tone	5 sec Pure tone
Rest	15 ms	15 ms	15 ms
test frequencies	50 Hz 63 Hz 80 Hz 100 Hz 125 Hz 160 Hz 200 Hz 250 Hz 315 Hz	400 Hz 500 Hz 630 Hz 800 Hz 1000 Hz 1250 Hz 1600 Hz 2000 Hz 2500 Hz	3150 Hz 4000 Hz 5000 Hz 6300 Hz 8000 Hz 10000 Hz 12500 Hz 16000 Hz 20000 Hz
start frequency	50 Hz	400 Hz	3150 Hz
straat test ear side	right	right	right
intensity	15 dB, 30 dB, 45 dB, 60 dB,75 dB, 90 dB, 105 dB	15 dB, 30 dB, 45 dB, 60 dB,75 dB, 90 dB, 105 dB	15 dB, 30 dB, 45 dB, 60 dB,75 dB, 90 dB, 105 dB
step size	15 dB	15 dB	15 dB


Lawrence Neal 5:24 PM

Hey all, FYI, YouTube URL (<https://www.youtube.com/watch?v=tXJj61EptSk>) for DAY 4 - APRIL 28 BCIs WITH fNIRS, STIMULATION & ERROR POTENTIALS - 19:00 | Wireless EEG and fNIRS recordings Patrick Reitner.. starts @ 9:51:42, I think it's worth watching again.

fNIRS Principal




Continuous Wave fNIRS Principal




$$\begin{bmatrix} \Delta[Hb] \\ \Delta[O_2Hb] \end{bmatrix} = (d)^{-1} \begin{bmatrix} \epsilon_{Hb, \lambda_1} & \epsilon_{O_2Hb, \lambda_1} \\ \epsilon_{Hb, \lambda_2} & \epsilon_{O_2Hb, \lambda_2} \end{bmatrix}^{-1} \begin{bmatrix} \frac{\Delta OD(\Delta t, \lambda_1)}{DPF(\lambda_1)} \\ \frac{\Delta OD(\Delta t, \lambda_2)}{DPF(\lambda_2)} \end{bmatrix}$$

HBO and HBR



- oxygenated hemoglobine (HBO)
- deoxygenated hemoglobine (HBR)
- SNR in general better for HBO
- clear visible pulsatile fluctuations (quality indication)


HBO Response




- averaged HBO responses from MI experiment:
 - dotted line: task onset
 - task duration: 10s
- initial dip (~2-3s)
- peak (~8s)
- undershoot

HBR Response

- averaged HBR responses from MI experiment:
 - dotted line: task onset
 - task duration: 10s




fNIRS and EEG Placement




g.Nautilus fNIRS: Hard Facts

- EEG: g.Nautilus
 - up to 64 channels
 - 24 bit
 - up to 500Hz sampling rate
- fNIRS: Artinis OctaMon
 - 8 channels
 - 50Hz sampling rate
 - 760nm and 850nm


MI Experiment



MI Experiment – Simulink Model




MI Experiment – HBO Response



MI Experiment – HBR Response


Left Hand Grasp MI

Right Hand Grasp MI



lower on the right hemisphere, this is what we call

MI Experiment - Results



So, and now I'm going to show you a slide, which



10:03:13 / 10:37:52




Recap

- g.Nautilus fNIRS:
 - up to 64 channel EEG + 8 channel fNIRS (2 wavelengths)
- HBO + HBR
- Hemodynamic response
- Neurovascular coupling
- robustness + accuracy improvement
- multitude of BCI applications


That's reason why we combine



10:04:46 / 10:37:52



The next slide/info from Day 4...



Optical Densities

Concentrations

Sample Counter

Signal Quality

NIRS Interface

Data Acquisition - Device Selection

SERVER LOCATION

IP Address:

127.0.0.1

Port:

50221

Find connected devices...

CONNECTED DEVICES

Device Name	Device Type	In Use
RD-100L09.02	gNautilus	no

Test

Up

Down

Last

Impedance Check

Calibrate

Select

Close

UB-0000.00.00

EEG Data

double

Data Type Conversion

g.Hlsys

g.tec medical engineering GmbH

Display


Constant



gNautilus


Ground

10:07:14 / 10:37:52

View diagnostics







Patrick Reintner

